

# EVAPORATING DISH (W024)

Category: Lab Tools



**Product Code:** W024

**Category:** LAB TOOLS

**Sub-Category:** GLASSWARE

- 
- An evaporating dish in engineering testing laboratories is used for evaporating solvents from a solution, leaving behind the solid residue. It is typically made of porcelain or heat-

resistant glass, designed to withstand high temperatures during the evaporation process, commonly used in chemical analyses or material testing.

---

## **ORDERING INFORMATION**

| Item                        | Code      |
|-----------------------------|-----------|
| EVAPORATING DISH - Ø 30 mm  | W024C015H |
| EVAPORATING DISH - Ø 70 mm  | W024C050H |
| EVAPORATING DISH - Ø 80 mm  | W024C075H |
| EVAPORATING DISH - Ø 100 mm | W025C140H |
| EVAPORATING DISH - Ø 125 mm | W025C300H |

---

*Keywords: evaporating dish, W023 evaporating dish, ALFA Testing Equipment, laboratory glassware, solvent evaporation tool, residue recovery dish, engineering lab equipment, high-temperature heating tool*